



# NB

6 10 15 20 25 30 35 40 45 50 55 60 65 70 75 80 85 90 95 100

- NUVE NB 5
- NUVE NB 9
- NUVE NB 15
- NUVE NB 20
- NUVE NB 25
- NUVE NB 30
- NUVE NB 35
- NUVE NB 40
- NUVE NB 45
- NUVE NB 50
- NUVE NB 55
- NUVE NB 60
- NUVE NB 65
- NUVE NB 70
- NUVE NB 75
- NUVE NB 80
- NUVE NB 85
- NUVE NB 90
- NUVE NB 95
- NUVE NB 100



# N-prize

|  | NB 9                                  | NB 20       |
|--|---------------------------------------|-------------|
| Tank Volume, litres                            | 9,5                                   | 21          |
| Useful Volume, litres                          | 7                                     | 15          |
| Temperature Range                              | Ambient Temp + 5°C / 99,9°C           |             |
| Temperature Sensor                             | Fe-Const                              |             |
| Control System                                 | N-Prime™ Programmable Microprocessor  |             |
| Displays                                       | Displays for Temperature and Time     |             |
| Temperature Set and Display Sensitivity        | 0,1°C                                 |             |
| Temperature variation @37°C                    | ± 0,1°C*                              |             |
| Temperature fluctuation @37°C                  | ± 0,01°C*                             |             |
| Programmable Alarm Limits                      | ± 0,5°C / 5°C                         |             |
| Timer  | 1 minute – 99,9 hours + Hold Position |             |
| Delayed Start Timer                            | 1 minute – 99,9 hours                 |             |
| Selectable Recording Frequency on Memory Stick | 1-999 minutes                         |             |
| Internal Material                              | Stainless Steel                       |             |
| External Material                              | Epoxy - Polyester Powder Coated Steel |             |
| Power Rating                                   | 800 W                                 | 1600 W      |
| Power Supply                                   | 230 V, 50/60 Hz                       |             |
| Internal Dimension (WxDxH) mm                  | 240x300x150                           | 505x300x150 |
| External Dimension (WxDxH) mm                  | 325x400x275                           | 590x405x275 |
| Packing Dimension (WxDxH) mm                   | 390x450x480                           | 650x450x480 |
| Net/Packed Weight Kg.                          | 9 / 12                                | 13 / 17     |

\*At 22°C ambient temperature and closed lid

## FACTORY FITTED OPTION

NB XXX WF Wi-Fi connection for N-Mobile™

## ACCESSORIES

|          |   |          |                                    |
|----------|---|----------|------------------------------------|
| K 04 286 | Plexiglass lid for NB 9 (Resists up to 60°C)  | A 08 051 | Tube rack 52 x Ø 13 mm, wide 70 mm |
| K 52 007 | Stainless steel lid for NB 9                  | A 08 050 | Tube rack 30 x Ø 16 mm, wide 70 mm |
| K 04 289 | Plexiglass lid for NB 20 (Resists up to 60°C) | A 08 021 | Tube rack 27 x Ø 18 mm, wide 70 mm |
| K 52 003 | Stainless steel lid for NB 20                 | A 08 049 | Tube rack 12 x Ø 30 mm, wide 70 mm |
| K 52 008 | Lid with 4 holes for NB 9 (Hole dia. 95 mm)   | R 01 113 | 2x1000ml Flask holder              |
| K 52 005 | Lid with 4 holes for NB 20 (Hole dia. 95 mm)  | R 01 116 | 4x500ml Flask holder               |
| K 52 006 | Lid with 6 holes for NB 20 (Hole dia. 95 mm)  |          |                                    |
| K 52 014 | Lid with 8 holes for NB 20 (Hole dia. 95 mm)  |          |                                    |



**NUVE SANAYİ MALZEMELERİ İMALAT VE TİCARET A.Ş.**

Saracalar Mah. Saracalar Kümeevleri No: 4/2 Akyurt 06750 ANKARA / TÜRKİYE t. +90 312 399 28 30 f. +90 312 399 21 97  
nuve.com.tr sales@nuve.com.tr

ISO 9001:2015  
ISO 13485:2016

